

ENgage with Every Moment

The PUGGLE™ Enteral Feeding Pump is designed to complement the lifestyles of people with enteral feeding requirements, so each person can live their best life.

From infants to seniors, the PUGGLE™ Pump's thoughtful, smart design reliably supports patients whether they are resting or on-the-go. Our versatile programs and continuous motor work together to enable better management of enteral feeding. Our system delivers the nutrition and energy that patients need – how they need it, and when they need it.



ENjoy Life Optimized

LCD Screen features a large color display for an at-a-glance view of key settings and progress

Battery provides 20 hours of power with a speedy 2.5 hour charging time

Intuitive Design uses simple buttons and screen instructions for easy-to-manage programming



Intermittent Mode allows programming for up to 8 feeds, then automatically starts the program

Motor features a quiet, continuous peristaltic rotary pump action

Precise Performance with a flow rate and volume accuracy rate of $\pm 5\%$ over one hour

ENgineered with Love

Peace-of-mind with thought-through, intuitive design and reliable performance.



Item Code	Description	Dimensions	Weight	Flow Rate Range	Volume Range	Accuracy*
PG001	PUGGLE™ Enteral Feeding Pump – AC Adapter and Pole Clamp Included	5.4"H x 5.0"W x 2.5"D	1.43 lbs.	1 to 400 mL/h in 1 mL increments	1 to 9,999 mL in 1 mL increments	± 5%

* Accuracy of feed at the specified rate over one hour.

Item Code	Description	Qty/Unit
PE90ST	PUGGLE™ Enteral Pump Spike Set, Anti-free Flow Valve, and ENFit® Connector and Transition Connector	30/Case
PE90B5T	PUGGLE™ Enteral Pump Set with 500 mL Bag, Anti-free Flow Valve, and ENFit® Connector and Transition Connector	30/Case
PE90B12T	PUGGLE™ Enteral Pump Set with 1200 mL Bag, Anti-free Flow Valve, and ENFit® Connector and Transition Connector	30/Case



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768 • Amsino.com

Phone: 800.632.6746 • Fax: 909.626.3888 • Email: info@amsino.com

Please refer to the Instructions for Use for additional information.

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. ENFit® is a registered trademark of GEDSA. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © Copyright 2021, Amsino International, Inc. All rights reserved.

EN-US-2021-4-30-DS