

ENgage with Every Moment

The PUGGLE™ Enteral Feeding Pump is designed to complement the lifestyles of people with enteral feeding requirements, so each person can live their best life.

From infants to seniors, the PUGGLE™ Pump's thoughtful, smart design reliably supports patients whether they are resting or on-the-go. Our versatile programs and continuous motor work together to enable better management of enteral feeding. Our system delivers the nutrition and energy that patients need – how they need it, and when they need it.

ENjoy Life Optimized

LCD Screen features a large color display for an at-a-glance view of key settings and progress

Battery provides 20 hours of power with a speedy 2.5 hour charging time

Intuitive Design uses simple buttons and screen instructions for easy-to-manage programming



Intermittent Mode allows programming for up to 8 feeds, then automatically starts the program

Motor features a quiet, continuous peristaltic rotary pump action

Precise Performance with a flow rate and volume accuracy rate of $\pm 5\%$ over one hour

ENgineered with Love

Peace-of-mind with thought-through, intuitive design and reliable performance.



Item Code	Description	Dimensions	Weight	Flow Rate Range	Volume Range	Accuracy*
PG001	PUGGLE™ Enteral Feeding Pump – AC Adapter and Pole Clamp Included	5.4"H x 5.0"W x 2.5"D	1.43 lbs.	1 to 400 mL/h in 1 mL increments	1 to 9,999 mL in 1 mL increments	± 5%

* Accuracy of feed at the specified rate over one hour.

Item Code	Description	Qty/Unit
PE90ST	PUGGLE™ Enteral Pump Spike Set, Anti-free Flow Valve, and ENFit® Connector and Transition Connector	30/Case
PE90B5T	PUGGLE™ Enteral Pump Set with 500 mL Bag, Anti-free Flow Valve, and ENFit® Connector and Transition Connector	30/Case
PE90B12T	PUGGLE™ Enteral Pump Set with 1200 mL Bag, Anti-free Flow Valve, and ENFit® Connector and Transition Connector	30/Case



Amsino International, Inc.
 708 Corporate Center Drive, Pomona, CA 91768 • Amsino.com
 Phone: 800.632.6746 • Fax: 909.626.3888
 Email: info@amsino.com

PUGGLE™ Instructional Video



Please refer to the Instructions for Use for additional information.

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. ENFit® is a registered trademark of GEDSA. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © Copyright 2021, Amsino International, Inc. All rights reserved.

EN-US-2021-4-30-DS