



## ***EN*joy Life Uninterrupted**

The PUGGLE™ Pump is designed to complement each EN patient's nutritional needs and lifestyle.



# ENjoy Life Optimized

The versatile PUGGLE™ Enteral Feeding Pump is designed to deliver nutritional solutions while maintaining the dignity and independent lifestyle for patients with enteral feeding needs. From infants to seniors, the PUGGLE Pump offers stationary and ambulatory applications. It's not only life-friendly, it's life-nourishing.

**Motor** – Features a quiet, continuous peristaltic rotary pump action

**LCD screen** – Large color display provides an at-a-glance view of key settings and progress

**Precise** – Delivers flow rate and volume feed at an accuracy rate of  $\pm 5\%$  over one hour



**Intermittent mode** – Allows programming for up to 8 feeds

**Battery** – Provides 20 hours of power with a speedy 2.5 hour charging time

**Intuitive** – Simple buttons and screen instructions make programming easy to manage



## ENgage with Every Moment

Peace-of-mind with thought-through, intuitive design and reliable performance.

### Nutrition & Life – Uninterrupted

The PUGGLE™ Pump's innovation works quietly in the background so patients can enjoy an optimized lifestyle.

### Thoughtful Functionality

The unique intermittent mode allows users to program up to 8 automatic feeds to deliver a seamless, yet secure, nutrition experience.

### All-Day Dependability

The 20-hour battery life and 2.5-hour charging time allows patients to focus on the NOW, instead of sitting out the day for a battery charge.

### Accommodating and Compact Design

The PUGGLE™ Pump enables patients to receive vital nutrition whether they're resting at home or on the go.

### Intuitive Technology

The PUGGLE™ Pump features a user-friendly interface that is secure and easy to manage.

### Ready, Set, Pump

Our quick and automatic priming system allows patients to get started in under 90 seconds so no time is wasted.

## ENgineered With Love

We believe that people with EN needs shouldn't have to fuss with a clunky device to receive their vital nutrients. They should be able to freely live an optimized life. That's why the PUGGLE™ Pump is engineered to give caregivers and patients peace of mind through an intuitive design and its reliable, results-driven performance.

For over 28 years, Amsino has been serving the needs of the healthcare industry in pursuit of improving patient care. Amsino develops, manufactures, and markets medical products that help improve the safety and effectiveness of patient care under a portfolio of recognized Amsino brands.



Quiet Motor



Robust Battery Life



Easy To Use



Flexible Settings



For Kids & Adults

To learn more visit [PUGGLEpump.com](https://PUGGLEpump.com)

The logo for AMSINO, featuring the word "AMSINO" in a bold, white, serif font, set against a dark blue rectangular background.

Advancing Healthcare Worldwide®

**Amsino International, Inc.**

708 Corporate Center Drive, Pomona, CA 91768 • [Amsino.com](http://Amsino.com)

**Phone:** 800.632.6746 • **Fax:** 909.626.3888 • **Email:** [info@amsino.com](mailto:info@amsino.com)

Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings, and precautions.

Product(s) described may not be licensed or available for sale in all countries. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © Copyright 2021, Amsino International, Inc. All rights reserved.

EN-US-2021-8-24-BRO